
Extra choices

For disposable electro-surgical pencil:

- Needle/ball/loop tip.
- 7, 8, 10, 12, 14, 16cm long electrode.

For disposable patient plate:

- 5m pre-attached cable.

*For other customization needs, please contact Hisern's sales for details.



HISERN

Web: www.hisern.com/en

E-mail: sales@hisern.com / info@hisern.com

Tel: +86 0755-86660662

© 2024 Zhejiang Haisheng Medical Device Co., Ltd. All rights are reserved.

Zhejiang Haisheng Medical Device Co., Ltd. reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

HISERNES00120240001



LinkedIn



Website


Electrosurgery Series

HISERN

Disposable Electrosurgical Pencil

Electrosurgical pencil is widely used in various surgical operations for electrocutting and electrocoagulation.


- Ergonomic design
- Lightweight
- Ease of use
- Button switch
- Improved safety



PN	Model	Type	Coating	Tip	Electrode	Cable	Pcs/Carton
H01009001	HEP4	Standard	Stainless Steel	Blade	4 cm	2.5m	100
H01009002	HEP4	Standard	Stainless Steel	Blade	4 cm	3m	100
H01009003	HEP4	Standard	Stainless Steel	Blade	4 cm	5m	100
H01009004	HEPS4	Standard	PTFE Coated	Blade	4 cm	2.5m	100
H01009005	HEPS4	Standard	PTFE Coated	Blade	4 cm	3m	100
H01009006	HEPS4	Standard	PTFE Coated	Blade	4 cm	5m	100
H01009007	HEPM	Smoke Evacuation	Stainless Steel	Blade	4 cm	2.5m	100
H01009008	HEPM	Smoke Evacuation	Stainless Steel	Blade	4 cm	3m	100
H01009009	HEPM	Smoke Evacuation	Stainless Steel	Blade	4 cm	5m	100
H01009010	HEPSM	Smoke Evacuation	PTFE Coated	Blade	4 cm	2.5m	100
H01009011	HEPSM	Smoke Evacuation	PTFE Coated	Blade	4 cm	3m	100
H01009012	HEPSM	Smoke Evacuation	PTFE Coated	Blade	4 cm	5m	100

Electrosurgical Pencil Tip Cleaner

PN	Size	Pcs/Carton
H00309001	4.8x4.8cm	400
H00309002	4.8x2.4cm	400



Patient Plate

Patient plate is used in electrosurgery to remove current from the patient safely


- Fills patient skin irregularities to improve the total electrical conduction area.
- Gel and adhesive foam are tolerated even by sensitive skin.
- Multiple configurations - available in pediatric and adult sizes.



Non-split / Monopolar



Split / Bipolar



Jack Connector



Non-REM Connector



REM Connector

PN	Model	Category	Split/Non-split	Cable	Connector	Pcs/Carton
H00509001	HPP122	Adult	Non-split/Monopolar	No cable	/	200
H00509002	HPP112	Adult	Non-split/Monopolar	3m	Jack	100
H00509003	HPP112	Adult	Non-split/Monopolar	3m	Non-REM	100
H00509004	HPP222	Adult	Split/Bipolar	No cable	/	200
H00509005	HPP212	Adult	Split/Bipolar	3m	Jack	100
H00509006	HPP212	Adult	Split/Bipolar	3m	Non-REM	100
H00509007	HPP212	Adult	Split/Bipolar	3m	REM	100
H00509008	HPP121	Pediatric	Non-split/Monopolar	No cable	/	200
H00509009	HPP111	Pediatric	Non-split/Monopolar	3m	Jack	100
H00509010	HPP111	Pediatric	Non-split/Monopolar	3m	Non-REM	100
H00509011	HPP221	Pediatric	Split/Bipolar	No cable	/	200
H00509012	HPP211	Pediatric	Split/Bipolar	3m	Jack	100
H00509013	HPP211	Pediatric	Split/Bipolar	3m	Non-REM	100
H00509014	HPP211	Pediatric	Split/Bipolar	3m	REM	100